

tacoTM Prep Bead Beater



Ease of Use

- Easy one-button operation

Efficient

- Homogenize various biological samples including animal tissues, plant tissues and microorganisms
- Process up to 16 samples per run

Flexible Application

- Adjustable cycle number for effective homogenization
- Homogenized products readily to be used for a wide range of downstream applications, such as nucleic acid extraction by **tacoTM** Nucleic Acid Extraction System

Specification	
Dimensions	30 mm × 25 mm × 24 mm, 11 kg
Power Supply	100~240 V, 50/60 Hz, AC
Sample Capacity	16
Cycle duration	13 sec/cycle or 40 sec/cycle
Sample Type	Biological samples (animal tissues, plant tissues or microorganisms)
Consumable	Preloaded tacoTM Prep tube - 100 pcs/package
	tacoTM Prep tube - 500 pcs/package
	tacoTM Prep steel beads - 3,000 pcs/package



GeneReach Biotechnology Corp

Tel: 886-4-2463-9869 E-mail: sales@genereach.com

No.19 Keyuan 2nd Rd. Central Taiwan Science Park, Taichung City, Taiwan